

## **SDS-PAGE Electrophoresis Buffer**

**Catalogue No.** AR0139

**Size** 1L

**Form** White powder (reconstitute with distilled water).

**Storage** Store at 4°C in dark for one year.

### **Protocol**

1. Pour a packet of SDS-PAGE powder (enough to prepare 1L Electrophoresis Buffer) into clean beaker.
2. Add distilled water until up to 900 ml, mix to dissolve. Keep the constant volume 1L and mix.
3. The remaining Electrophoresis buffer can be stored at room temperature for two weeks.