

### Rat Fibrinogen Total Antigen Assay ELISA Kit

Catalog No.	CSI20063A CSI20063B	1 x 96 tests 5 x 96 tests
Intended Use:	to insoluble fibrin by thrombin in the Hepatic expression of fibrinogen inc	that circulates in the blood and is converted final step of the coagulation cascade. reases two to four hundred fold during the acute phase n. Elevated fibrinogen levels are correlated with clerosis.
Assay Principle:	plate. After appropriate washing ste antibody binds to the captured prote polyclonal antibody is reacted with a substrate is used for color developm	purified capture antibody coated on the microtiter os, biotin labeled polyclonal anti-rat fibrinogen primary in. Excess antibody is washed away and bound vidin conjugated to horseradish peroxidase. TMB tent at 450nm. A standard calibration curve is prepared red using dilutions of rat fibrinogen. Color development of fibrinogen in the samples.
Reagents Provided:	<ul> <li>\$96-well microtiter strip plate (8X12 Fibrinogen capture antibody coated</li> <li>\$5X Diluent:</li> <li>1 bottle of 50ml; bring to 1X using D</li> <li>\$10X Wash Buffer:</li> <li>1 bottle of 50ml; bring to 1X using D</li> <li>\$Rat fibrinogen antigen standard:</li> <li>1 vial of lyophilized standard</li> <li>\$Anti-rat fibrinogen primary antibody</li> <li>1 vial of lyophilized biotin labeled po</li> <li>\$Avidin peroxidase secondary reage</li> <li>1 vial of concentrated HRP labeled a</li> <li>\$TMB substrate solution:</li> <li>1 bottle of 10ml solution</li> <li>\$Stop solution:</li> <li>1 bottle of 6ml 1M sulfuric acid</li> </ul>	l water I water r: lyclonal antibody ent:
Storage and Stability:	microtiter strips in the pouch with de stored at -70° C for later use. <b>DO No</b>	t 4°C. Store unopened plate and any unused siccant. Reconstituted standards and primary may be <b>DT</b> freeze/thaw the standards and primary antibody more bonents must be stored at 4°C. Kit should be used no



**Cell Sciences<sup>®</sup>** Neponset Valley Tech Park 480 Neponset St., Bldg. 12A Canton, MA 02021 Toll Free:888 769-1246Phone:781 828-0610Fax:781 828-0542



Reagents and Equipment Required:	<ul> <li>1-channel pipettes covering 0-10µl and 200-1000µl</li> <li>12-channel pipette covering 30-300µl</li> <li>Paper towels or kimwipes</li> <li>50ml tubes, 1.5ml centrifuge tubes</li> <li>1N H<sub>2</sub>SO<sub>4</sub></li> <li>DI water</li> <li>Magnetic stirrer and stir-bars</li> <li>Plastic containers with lids</li> <li>Microtiter plate spectrophotometer operable at 450nm</li> <li>Microtiter plate shaker with uniform horizontally circular movement up to 300rpm. (OPTIONAL)</li> </ul>
Warnings:	Warning – Avoid skin and eye contact when using TMB One substrate solution since it may be irritating to eyes, skin, and respiratory system. Wear safety goggles and gloves.
Precautions:	<ul> <li>•DO NOT mix any reagents or components of this kit with any reagents or components of any other kit. This kit is designed to work properly as provided.</li> <li>•DO NOT pipette reagents by mouth.</li> <li>•Always pour substrate out of the bottle into a clean test tube. DO NOT pipette out of the bottle as you could contaminate the substrate.</li> <li>•Keep plate covered except when adding reagents, washing, or reading.</li> <li>•All kit components must be kept refrigerated (4°C).</li> <li>•DO NOT smoke, drink, or eat in areas where specimens or reagents are being handled.</li> </ul>
Preparation of Reagents:	<ul> <li>•Diluent concentrate: The diluent supplied in a 5X concentrate and must be diluted 1:5 with deionized water for use with the kit.</li> <li>•Wash buffer concentrate: The wash buffer supplied in a 10X concentrate and must be diluted 1:10 with deionized water for use with the kit.</li> </ul>
Specimen Collection:	The assay measures total rat fibrinogen in the 3.125-800 ng/ml range. Samples giving rat fibrinogen levels above 800ng/ml should be diluted in 1X diluent before use. Normal plasma samples need to be diluted between 1:10,000 and 1:50,000 in 1X diluent for the values to be within linear range of the standard curve.
Assay Procedure:	Perform assay at room temperature. Vigorously shake plate (300rpm) at each step of the assay. If a microtiter plate shaker is not available then keep the plate on a flat surface for 60 minutes at each step instead of 30 minutes. Note: when the assay is performed without shaking the plate, the final absorbance values at 450nm will be lower than when the assay is performed using a plate shaker.



Toll Free:888 769-1246Phone:781 828-0610Fax:781 828-0542



Preparation of Standard:

Reconstitute standard vial with 5 ml of 1X diluent to give a 800ng/ml solution. Dilution table for preparation of rat fibrinogen standards:

Fibrinogen	Dilutions
concentration	
(ng/ml)	
800	Straight from the vial
000	ou aight norm and via
400	500µl (1X Diluent) +
	500µl (800ng/ml)
200	500µl (1X Diluent) +
	500µl (400ng/ml)
100	500µl (1X Diluent) +
	500µl (200ng/ml)
50	500µl (1X Diluent) +
	500µl (100ng/ml)
25	500µl (1X Diluent) +
	500µl (50ng/ml)
12.5	500µl (1X Diluent) +
	500µl (25ng/ml)
6.25	500µl (1X Diluent) +
	500µl (12.5ng/ml)
3.125	500µl (1X Diluent) +
	500µl (6.25ng/ml)
0	500µl (1X Diluent)
	Zero point to determine
	background

## NOTE: DILUTIONS FOR THE STANDARD CURVE AND ZERO STANDARD MUST BE MADE AND APPLIED TO THE PLATE IMMEDIATELY.

### Standard and Unknown Addition:

Remove microtiter plate from bag. Add 100µl standards in duplicate and unknowns to wells. Carefully record position of standards and unknowns. Shake plate at 300rpm for 30 minutes. Wash wells three times with 300µl wash buffer. Remove excess wash by gently tapping plate on paper towel or kimwipe.

### Primary Antibody Addition:

Add 10ml of 1X diluent directly to the primary antibody vial and agitate gently to completely dissolve contents. Add 100µl to all wells. Shake plate at 300rpm for 30 minutes. Wash wells three times with 300µl wash buffer. Remove excess wash by gently tapping plate on paper towel or kinwipe.

### Secondary Antibody Addition:

Dilute 2.5µl into 10ml of 1X diluent and mix well and add 100µl to all wells. Shake plate at 300rpm for 30 minutes. Wash wells three times with 300µl wash buffer. Remove excess wash by gently tapping plate on paper towel or kinwipe.



**Cell Sciences**<sup>®</sup> Neponset Valley Tech Park 480 Neponset St., Bldg. 12A Canton, MA 02021 Toll Free:888 769-1246Phone:781 828-0610Fax:781 828-0542

# cellsciences.com

### Substrate Incubation:


Add 100 $\mu$ I TMB substrate to all wells and shake plate for 10-20 minutes. Substrate will change from colorless to different strengths of blue. Quench reaction by adding 50 $\mu$ I of 1N H<sub>2</sub>SO<sub>4</sub> stop solution to all wells when samples are visually in the same range as the standards. Add stop solution to wells in the same order as substrate upon which color will change from blue to yellow. Mix thoroughly and read final absorbance values at 450nm. For best results read plate immediately.

### Measurement:

Set the absorbance at 450nm in a microtiter plate spectrophotometer. Measure the absorbance in all wells at 450nm,  $A_{450}$ . Subtract zero point from all standards and unknowns to determine corrected absorbance ( $A_{450}$ ).

### Assay Calibration:

Plot  $A_{450}$  against the amount of rat fibrinogen in the standards. Fit a straight line through the points using a linear fit procedure. The amount of total rat fibrinogen in the unknowns can be determined from this curve.





**Cell Sciences<sup>®</sup>** Neponset Valley Tech Park 480 Neponset St., Bldg. 12A Canton, MA 02021 Toll Free:888 769-1246Phone:781 828-0610Fax:781 828-0542



ExpectedFibrinogen is present in normal rat plasma at a concentration of 3.1 mg/ml and varies by age and<br/>diet. Normal plasma samples need to be diluted between 1:10,000 and 1:50,000 in<br/>1XDiluent for the values to be within linear range of the standard curve

**Disclaimer:** This information is believed to be correct but does not claim to be all-inclusive and shall be used only as a guide. The supplier of this kit shall not be held liable for any damage resulting from handling or from contact with the above product.

#### Example of the plate layout:

	1	2	3	4	5	6	7	8	9	10	11	12
		3.125	6.25	12.5	25	50	100	200	400	800		
Α	0	ng/ml										
		3.125	6.25	12.5	25	50	100	200	400	800		
в	0	ng/ml										
с												
D												
Е												
F												
G												
н												
96 N	96 well plate											
	Otendarda, Otenalla											

Standards: 20 wells Samples: 76 wells

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



**Cell Sciences**<sup>®</sup> Neponset Valley Tech Park 480 Neponset St., Bldg. 12A Canton, MA 02021 Toll Free:888 769-1246Phone:781 828-0610Fax:781 828-0542