

Human BFAR cDNA Clone



Sino Biological Inc.
Biological Solution Specialist

Catalog Number: HG11334-M

General Information

Gene Bank Ref. ID : NM_016561.2

cDNA Size: 1353

Gene Name Synonym :

BAR; RNF47; BFAR

Source: Human

Description

cDNA Description :

ORF Clone of Homo sapiens bifunctional apoptosis regulator DNA.

Vector :

Commercial pMD vector. (This vector is produced by TAKARA.

More information: <http://www.takara.com.cn/explain/L/D103A.pdf>)

Sequence Description :

Identical with the Gene Bank Ref. ID sequence.

Shipping

1. 5 - 10µg Plasmid DNA is shipped on Whatman FTA elute card (Catalog Number: WB120410) at room temperature.
2. It would take 3 - 5 days for delivery typically. For bulk orders, additional 1 - 2 weeks would be needed.

Usage Protocol

1. Remove the marked edge in the center of the spot.
2. Transfer the filter paper to a 1.5ml microfuge tube containing 50-100 µl of sterile H₂O. Ensure the filter is completely immersed in the H₂O by briefly centrifuging the tube for 10 seconds.
3. Transfer the tube to a heating block at 95°C for 15 - 30 min.
4. At the end of the incubation period remove the sample from the block and pulse vortex, or gently tap, the sample approximately 60 times.
5. Briefly centrifuge for 30 seconds, to separate the matrix from the elute. The elute now contains the purified DNA.

cDNA Sequence and Quality Control

Sequence :

```
ATGGAGGAACCTCAGAAAAGCTATGTGAACACAATGGACCTTGAGAGAGATGAA
CCTCTCAAAGCACCCGGCCCTCAGATTTCTGTTAGTGAATTTCTTGCCACTGCT
GCTACGACATCCTGGTTAACCCACCACCTTGAAGTGTGGGCACAGCTTCTGCC
GTCAGTGCCTTGTCTTATGGTGGGCATCTTCAAAGAAAACAGAAATGTCAGAA
GCAGAGAAAAATGGGAAGGTTTCCCCAAAGTCAGATTCTCCTCAGGGATGCCA
TTGAAAAGTTATTTCTGATGCCATTAGACTGAGATTGAAGACATTCAGCAGAA
TAATGACATAGTCCAAAGTCTTGCAGCCTTTCAGAAATATGGGAATGATCAGAT
CCTTAGCTCCTAACACAGGCCGAGCGAATCAGCAGATGGGAGGGGGATTCTT
TTCCGGTGTGCTCACAGCTTTAACTGGAGTGGCAGTGGTCTGCTCGTCTATCA
CTGGAGCAGCAGGGAATCTGAACACGACCTCCTGGTCCACAAGGCTGTGGCCA
AATGGACGGCGGAAGAAGTTGCTCTGGCTGGAGCAGCTGGGCCCTTGGGC
ATCTTTTACAGGGAAGGTTTTATCTGAACGAGTAAATGGAAGGTTGCTTTTA
ACTTTGACAGAGGAAGAATTTTCCAAGACGCCCTATACCATAGAAAACAGCAGC
CACAGGAGAGCCATCCTCATGGAGCTAGAACGTGTCAAAGCATTAGGCGTGAA
GCCCCCAGAAATCTCTGGGAATATAAGGCTGTGAACCCAGGCAGGTCCTGT
TCCTGCTATACGCCCTCAAGAGCTCCCCCAGGCTGAGTCTGCTCTACCTGTACC
TGTTTGACTACACCGACACCTTCTACCTTTCATCCACACCATCTGCCCTCTGCA
AGAAGACAGCTCTGGGGAGGACATCGTCACCAAGCTTCTGGATCTTAAGGAGC
CTACGTGGAAGCAGTGGAGAGAGTTCTGTTCAATACTCCTTCCATCACC
AGCTGATTGCTGAGTTTGGTGGGACTGGTTGGAGGTCATTACTGGACATCAC
GGTTTCTCATCATCAATGCTATGTTACTCTCAGTTTCTGGAATTATTCTCCTTTGG
AGAATCTGGTCGAGAAGTGAAGTGAAGACCGTGCCTCAGAGGATGTGGAGCCA
TTTCTGGAAAGTATCAACGCAGGGGCTTTTTGTGGCCATGTTCTGGCCCTCAT
CCCTCAGTTTGTGCAACTGTTTGTGTTTACTGGGCCCTGTACTTTAACCAATTA
TTAACATTGATCTTGTGGTCAAGGAACCTCCGGCGGCTGGAACCCAGGTGTTGT
GA
```

Quality Control:

The sequencing result is available upon request.