

Home sapiens PCBP1 cDNA Clone



Sino Biological Inc.
Biological Solution Specialist

Catalog Number: HG11628-M

General Information

Gene Bank Ref. ID : NM_006196.3

cDNA Size: 1071

Gene Name Synonym :

HNRPX; HNRPE1; hnRNP-X; hnRNP-E1; PCBP1

Source: Human

Description

cDNA Description :

ORF Clone of Homo sapiens poly(rC) binding protein 1 DNA.

Vector :

Commercial pMD vector.(Produced by TAKARA. More information:

<http://www.takara.com.cn/explain/L/pMD18-T Simple.pdf>)

Sequence Description :

Identical with the Gene Bank Ref. ID sequence.

Shipping

1. 5 - 10µg Plasmid DNA is shipped on Whatman FTA elute card (Catalog Number: WB120410) at room temperature.
2. It would take 3 - 5 days for delivery typically. For bulk orders, additional 1 - 2 weeks would be needed.

Usage Protocol

1. Remove the marked edge in the center of the spot.
2. Transfer the filter paper to a 1.5ml microfuge tube containing 50-100 µl of sterile H₂O. Ensure the filter is completely immersed in the H₂O by briefly centrifuging the tube for 10 seconds.
3. Transfer the tube to a heating block at 95°C for 15 - 30 min.
4. At the end of the incubation period remove the sample from the block and pulse vortex, or gently tap, the sample approximately 60 times.
5. Briefly centrifuge for 30 seconds, to separate the matrix from the elute. The elute now contains the purified DNA.

cDNA Sequence and Quality Control

Sequence :

```
ATGGATGCCGGTGTGACTGAAAGTGGACTAAATGTGACTCTCACCATTGGCTT
CTTATGCACGGAAAGGAAGTAGGAAGCATCATTGGGAAGAAAGGGAGTCGGT
TAAGAGGATCCGCGAGGAGAGTGGCGCGGGATCAACATCTCGGAGGGGAAT
TGTCCGGAGAGAATCATCACTCTGACCGGCCCAACCAATGCCATCTTTAAGGCT
TTCCGCTATGATCATCGACAAGCTGGAGGAAGATATCAACAGCTCCATGACCAAC
AGTACCGCGGGCCAGCAGGCCCGGTACCCTGAGGCTGGTGGTGCCGGCCA
CCCAGTGC GGCTCCCTGATTGGGAAAGCGGGTGT AAGATCAAAGAGATCCGC
GAGAGTACGGGGGCGCAGGTCCAGGTGGCGGGGGATATGCTGCCAACTCCA
CCGAGCGGGCCATCACCATCGCTGGCGTGCCGCACTGTCCACCGAGTGTGTC
AAGCAGATTTGCCTGGTCATGCTGGAGACGCTCTCCAGTCTCCGCAAGGGAG
AGTCATGACCATTCCGTACCAGCCATGCCGGCCAGCTCCCCAGTCATCTGCG
CGGGCGGCCAAGATCGGTGCAGCGACGCTGCCGGCTACCCCATGCCACCA
TGACCTGGAGGACCACCTCTAGATGCCTACTGATTAAGGACAACACCCAT
TTCTCCGCTCGATCTGGCCAAGCTGAACCAGGTGGCAAGACAACAGTCTCACTT
TGCCATGATGCACGGCGGGACCGGATTCGCCGGAATTGACTCCAGCTCTCCAG
AGGTGAAAGGCTATTGGGCAAGTTTGGATGCATCTACTCAAACCACCCATGAAC
TCACCATCCAAATAACTTAATTGGCTGCATAATCGGGCGCCAAGCGCCAACA
TTAATGAGATCCGCCAGATGTCGGGGGCCAGATCAAATGGCCAACCCAGTG
GAAGGCTCCTCTGGTAGGCAGGTTACTATCACTGGCTCTGCTGCCAGTATTAGT
CTGGCCAGTATCTAATCAATGCCAGGCTTTCTCTGAGAAGGGCATGGGGTG
CAGCTAG
```

Quality Control:

The sequencing result is available upon request.