Homo sapiens EPYC cDNA Clone

Catalog Number: HG13855-G-F



General Information

Gene: epiphycan

Official Symbol: **EPYC**

Synonym: DSPG3, PGLB, Pg-Lb, SLRR3B, EPYC

Source: Homo sapiens

cDNA Size: 969

RefSeq: BC030958

Description

Please refer to the label on the tube Lot:

Sequence Description:

Identical with the Gene Bank Ref. ID sequence.

Restriction site: Kpnl + Xhol

Vector:

pCMV / hygro - FLAG

Shipping carrier:

Each tube contains approximately 10 µg of lyophilized plasmid.

Storage:

The lyophilized plasmid can be stored at ambient temperature for three months.

Quality control:

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

Sequencing primer list:

pcDNA3-L: 5' CTAGAGAACCCACTGCTTACTGGC 3' pcDNA3-R: 5' GGCAACTAGAAGGCACAGTCGAGG 3' Or

Forward T7: 5' TAATACGACTCACTATAGGG 3' ReverseBGH: 5' TAGAAGGCACAGTCGAGG 3'

pcDNA3-L and pcDNA3-R are designed by Sino Biological Inc. Customers can order the primer pair from any oligonucleotide supplier.

Plasmid Resuspension protocol

- 1. Centrifuge at $5,000 \times g$ for 5 min.
- 2. Carefully open the tube and add 100 µl of sterile water to dissolve the DNA.
- 3.Close the tube and incubate for 10 minutes at room temperature.
- 4.Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than $5000\times q$.
- 5. Store the plasmid at $-20 \,^{\circ}$ C.

The plasmid is ready for:

- Restriction enzyme digestion
- PCR amplification
- E. coli transformation
- DNA sequencing

E.coli strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5α and TOP10F'.

Fax:+86-10-51029969

● Tel:+86- 400-890-9989 ● http://www.sinobiological.com

Homo sapiens EPYC cDNA Clone

Catalog Number: HG13855-G-F




Vector Information

All of the pCMV vectors are designed for high-level stable and transient expression in mammalian hosts. High-level stable and non-replicative transient expression can be carried out in most mammalian cells. The vectors contain the following elements:

- · Human cytomegalovirus immediate-early (CMV) promoter for high-level expression in a wide range of mammalian cells.
- Hygromycin resistance gene for selection of mammalian cell
- A Kozak consensus sequence to enhance mammalian expression.

pCMV / hygro - FLAG
5681bp
Mammalian Expression Vector
Constitutive, Stable / Transient
CMV
Ampicillin
Hygromycin
FLAG

Physical Map of Plasmid:



http://www.sinobiological.com/Mammalian-Expression-Vectors-a-4666.html for the vector sequence.

Fax:+86-10-51029969