Mus musculus IL3Ra cDNA Clone

Catalog Number: MG50810-G



General Information

Gene: interleukin 3 receptor, alpha chain

Official Symbol: IL3RA

Synonym: CD123, SUT-1, CDw123

Source: Mus musculus

cDNA Size: 1191

RefSeq: NM_008369.1

Description

Lot: Please refer to the label on the tube

Sequence Description:

Identical with the Gene Bank Ref. ID sequence except for four point mutations: 338 A/G resulting in the amino acid asp substitution by gly, 367 C/A resulting in the amino acid leu substitution by met, 882 T/C and 1020 G/T not causing the amino acid variation.

Vector:

pGEM-T

Shipping carrier:

Each tube contains approximately 10 µg of lyophilized plasmid.

Storage:

The lyophilized plasmid can be stored at ambient temperature for three months.

Quality control:

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

Sequencing primer list:

M13-47: 5' GCCAGGGTTTTCCCAGTCACGAC 3'

RV-M: 5' GAGCGGATAACAATTTCACACAGG 3'

Other M13 primers can also be used as sequencing primers.

Plasmid Resuspension protocol

- 1. Centrifuge the tube for 5~10 min at 4,000 rpm.
- 2.Carefully open the tube and add 100 µl of sterile water to dissolve the DNA.
- 3. Close the tube and incubate for 10 minutes at room temperature.
- 4. Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than 4000 rpm.
- 5.Store the plasmid at -20 °C.

The plasmid is ready for:

- · Restriction enzyme digestion
- PCR amplification
- · E. coli transformation
- DNA sequencing

E.coli strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. DH5 α, TOP10, JM109.

Mus musculus IL3Ra cDNA Clone


Catalog Number: MG50810-G



Vector Information

The pGEM-T vector is a high-efficiency TA cloning vector which contains multiple cloning sites as shown below. The pGEM-T vector is 3.0kb in size and contains the amplicin resistance gene for selection. The coding sequence was inserted by TA cloning.

Physical Map of pGEM-T:



* Please refer to http://www.sinobiological.com/Vector-pGEM-T-a-1636.html for the vector sequence.