




Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK Telephone: +44 (0) 1223 755950 - Fax: +44 (0) 1223 755951 - E-Mail: info@abbexa.com

2K Plus DNA Marker (8 linear double-stranded DNA bands)

Catalogue No.:abx098046



2K Plus DNA Marker is composed of eight linear double-stranded DNA bands. This ready-to-use product is premixed with 1x DNA loading buffer. Based on experiment requirement, 5 µl marker can be directly loaded for gel electrophoresis. It provides simple operation and clear image.

Target: 2K Plus DNA Marker (8 linear double-stranded DNA bands)

Storage: Store at -20°C. Stable for 12 months from date of receipt.

Buffer: 10 mM Tris-HCl (pH 8.4), 10 mM EDTA, 0.02% bromophenol blue, 5% glycerol.

Note: This product is for research use only.