


DATASHEET

Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK Telephone: +44 (0) 1223 755950 - Fax: +44 (0) 1223 755951 - E-Mail: info@abbexa.com

12-120 kDa Protein Marker

Catalogue No.:abx098114



12-120 kDa Protein Marker. Protein Marker II is composed of eight unstained recombinant proteins 20 kDa, 30 kDa, 40 kDa, 50 kDa, 60 kDa, 80 kDa, 100 kDa, 120 kDa and one blue prestained recombinant protein 12 kDa. Prestained protein is obvious to monitor electrophoresis. 50 kDa band is a double intensity band for convenient identification of protein molecular weight. MW range from 12 kDa to 120 kDa. Ready-to-use format, no heating and reductant are required.

Target: 12-120 kDa Protein Marker

Storage: Store at -20°C. Stable for 12 months from date of receipt.

Buffer: 50 mM Tris-HCl pH 6.8, 5 mM EDTA, 10 mM DTT, 10% glycerol, 1% SDS, 0.01% phenol red.

Note: This product is for research use only.