

Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK Telephone: +44 (0) 1223 755950 - Fax: +44 (0) 1223 755951 - E-Mail: info@abbexa.com

DATASHEET

Rat Heat Shock 70 kDa Protein 1A (HSPA1A) ELISA Kit

Catalogue No.:abx155650



Rat HSPA1A ELISA Kit is a sandwich ELISA kit for use with Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids. This assay has high sensitivity and excellent specificity for detection of Heat Shock 70 kDa Protein 1A (HSPA1A)

No significant cross-reactivity or interference between Heat Shock 70 kDa Protein 1A (HSPA1A) and analogues was observed.

Please note that this kit is also available as a CLIA Kit <u>abx492397</u>.

Target:	Heat Shock 70 kDa Protein 1A (HSPA1A)
Reactivity:	Rat
Tested Applications:	ELISA
Recommended dilutions	Optimal dilutions/concentrations should be determined by the end user.
Immunogen:	Ala2-Glu386
Test Range:	1.56 ng/ml - 100 ng/ml
Sensitivity:	< 0.59 ng/ml
Validity:	The validity for this kit is 6 months.
Storage:	Store at 2°C to 8°C upon receipt.
Stability:	The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5% within the expiration date under appropriate storage conditions. To minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout.
Swiss Prot:	PODMWO

abbexa 🍊

DATASHEET

Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK Telephone: +44 (0) 1223 755950 - Fax: +44 (0) 1223 755951 - E-Mail: info@abbexa.com

GenelD:	<u>108348108</u>
Gene Symbol:	HSPA1A
Ensembl:	ENSRNOG0000045654
Standard Form:	Lyophilized
ELISA Detection:	Colorimetric
ELISA Type:	Sandwich
ELISA Data:	Quantitative
Sample Type:	Serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Note:	This product is for research use only.
	The range and sensitivity is subject to change. Please contact us for the latest product information. For accurate results, sample concentrations must be diluted to mid-range of the kit. If you require a
	specific range, please contact us in advance or write your request in your order comments.
	Please note that our ELISA and CLIA kits are optimised for detection of native samples, rather than recombinant proteins. We are unable to guarantee detection of recombinant proteins, as they may
	have different sequences or tertiary structures to the native protein.