

DATASHEET

Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK Telephone: +44 (0) 1223 755950 - Fax: +44 (0) 1223 755951 - E-Mail: info@abbexa.com

Inter Alpha-Globulin Inhibitor H4 (ITIH4) Antibody Pair

Catalogue No.:abx370764

Inter Alpha-Globulin Inhibitor H4 (ITIH4) Antibody Pair for use in Sandwich ELISA assay development. This antibody pair contains:

- Inter Alpha-Globulin Inhibitor H4 (ITIH4) capture antibody,
- Inter Alpha-Globulin Inhibitor H4 (ITIH4) biotin-conjugated detection antibody,
- Inter Alpha-Globulin Inhibitor H4 (ITIH4) standard.

It is recommended to use this antibody pair with abx098958 Antibody Pair Support Kit (Sandwich Method).

| Target: | Inter Alpha-Globulin Inhibitor H4 (ITIH4) |
|---------------------------------|--|
| Reactivity: | Rat |
| Tested Applications: | ELISA |
| Recommended dilutions: | Dilute the Capture Antibody with Coating Buffer. Dilute the biotin-conjugated Detection Antibody with Detection Antibody Diluent. Optimal dilutions/concentrations should be determined by the end user. |
| Form: | Liquid (Capture Antibody and Detection Antibody) |
| Reconstitution: | Reconstitute the standard with Standard Diluent. The volume, and therefore standard concentration, should be determined by the end user. |
| Storage: | Store at 2 to 8 °C for up to one month. Aliquot and store at -80 °C for up to one year. Avoid repeated freeze/thaw cycles. All solutions should be made fresh before the experiment. |
| Standard Form: | Lyophilized |
| ELISA Type: | Sandwich |
| Capture Antibody Conjugation: | Unconjugated |
| Detection Antibody Conjugation: | Biotin |
| Buffer: | The capture and detection antibody both contain 0.1% sodium azide. |



DATASHEET

Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK Telephone: +44 (0) 1223 755950 - Fax: +44 (0) 1223 755951 - E-Mail: info@abbexa.com

Directions for use:

Bring all components to room temperature (18-25°C) and briefly spin or centrifuge the vials before use. Working solutions should be prepared and used immediately.

Recommended Procedure:

- 1. Dilute the Capture Antibody to working concentration using Coating Buffer. Immediately coat the 96-well plate with diluted Capture Antibody (100 μ l per well). Seal the plate and incubate at 4 °C overnight or at 37 °C for 2 hours
- 2. Aspirate the wells and wash with Wash Buffer (350 µl per well) and allow to soak for 1-2 min. Remove the liquid by inverting and tapping the plate on to absorbent paper.
- 3. Block the plate with Blocking Buffer (200 µl per well) at 37 °C for 1.5 hours.
- 4. Repeat the aspiration/wash process in Step 2.
- 5. Add 100 μ l of standards or sample into the appropriate wells. Cover with a plate sealer and incubate at 37 °C for 1 hour.
- 6. Repeat the aspiration/wash process in Step 2.
- 7. Add appropriately diluted biotin-conjugated Detection Antibody (100 μ l per well). Cover the plate with a new plate sealer and incubate at 37 °C for 1 hour.
- 8. Repeat the aspiration/wash process in Step 2.
- 9. Add appropriately diluted Streptavidin HRP (100 μ l per well). Cover the plate with a new plate sealer and incubate at 37 °C for 30 min.
- 10. Repeat the aspiration/wash process in Step 2.
- 11. Add Substrate Solution (90 µl per well). Cover the plate with a new plate sealer and incubate at 37 °C for 10-20 min. Keep the plate in the dark and avoid exposure to light.
- 12. Add Stop Solution (50 µl per well). Tap the side of the plate to ensure thorough mixing.
- 13. Measure the absorbance immediately using a microplate reader set at 450 nm.

This product is for research use only.

Note: