

DATASHEET

Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK Telephone: +44 (0) 1223 755950 - Fax: +44 (0) 1223 755951 - E-Mail: info@abbexa.com

SLC26A8 siRNA

Catalogue No.:abx933786

siRNA to inhibit SLC26A8 expression using RNA interference.

This product is provided as two 5 nmol vials (10 nmol), three 5 nmol vials (15 nmol) or 2x three 5 nmol vials (30 nmol) of lyophilized siRNA oligo duplexes. Each vial contains slightly different sequences to ensure full knockout of the gene. The duplexes can be transfected individually or pooled together to achieve knockdown of the target gene, which is most commonly assessed by qPCR or western blot. The siRNA oligos are also chemically modified (2'-OMe) for increased stability and enhanced knockdown in vitro and in vivo.

Target: SLC26A8

Reactivity: Mouse

Host: Synthetic

Tested Applications: RNAi

Purity: > 97%

Form: Lyophilized

Specificity: SLC26A8 siRNA (Mouse) is a target-specific 19-23 nt siRNA oligo duplexes designed to knock down

gene expression.

Storage: Shipped at 4 °C. Store at -20 °C for up to one year.

Swiss Prot: Q8R0C3

GenelD: <u>224661</u>

Gene Symbol: SLC26A8

Directions for use: 1. Transfect with 100 nM siRNA 48 to 72 hours prior to cell lysis.

2. Before resuspending, briefly centrifuge the tube to ensure the lyophilized siRNA is at the bottom of

the tube.

3. Resuspend the siRNA oligos to an appropriate concentration with DEPC water. Each vial is suitable for 250 transfections in a 24 well plate (20 pmol for each well).

Quality Control: Oligonucleotide synthesis is monitored base by base through trityl analysis to ensure appropriate

coupling efficiency. The oligo is subsequently purified by affinity-solid phase extraction. The annealed RNA duplex is further analyzed by mass spectrometry to verify the exact composition of the duplex.

Each lot is compared to the previous lot by mass spectrometry to ensure maximum lot-to-lot

consistency.



DATASHEET

Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK Telephone: +44 (0) 1223 755950 - Fax: +44 (0) 1223 755951 - E-Mail: info@abbexa.com

Note: This product is for research use only.