## abbexa 👨

## **DATASHEET**

Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK Telephone: +44 (0) 1223 755950 - Fax: +44 (0) 1223 755951 - E-Mail: info@abbexa.com

## **ZSCAN4C siRNA**

Catalogue No.:abx941228

siRNA to inhibit ZSCAN4C expression using RNA interference.

This product is provided as two 5 nmol vials (10 nmol), three 5 nmol vials (15 nmol) or 2x three 5 nmol vials (30 nmol) of lyophilized siRNA oligo duplexes. Each vial contains slightly different sequences to ensure full knockout of the gene. The duplexes can be transfected individually or pooled together to achieve knockdown of the target gene, which is most commonly assessed by qPCR or western blot. The siRNA oligos are also chemically modified (2'-OMe) for increased stability and enhanced knockdown in vitro and in vivo.

Target: ZSCAN4C

Reactivity: Mouse

Host: Synthetic

Tested Applications: RNAi

**Purity:** > 97%

Form: Lyophilized

Specificity: ZSCAN4C siRNA (Mouse) is a target-specific 19-23 nt siRNA oligo duplexes designed to knock down

gene expression.

**Storage:** Shipped at 4 °C. Store at -20 °C for up to one year.

Swiss Prot: Q80VJ6

GenelD: <u>245109</u>

Gene Symbol: ZSCAN4C

Directions for use: 1. Transfect with 100 nM siRNA 48 to 72 hours prior to cell lysis.

2. Before resuspending, briefly centrifuge the tube to ensure the lyophilized siRNA is at the bottom of

the tube.

3. Resuspend the siRNA oligos to an appropriate concentration with DEPC water. Each vial is suitable for 250 transfections in a 24 well plate (20 pmol for each well).

**Quality Control:** Oligonucleotide synthesis is monitored base by base through trityl analysis to ensure appropriate

coupling efficiency. The oligo is subsequently purified by affinity-solid phase extraction. The annealed RNA duplex is further analyzed by mass spectrometry to verify the exact composition of the duplex.

Each lot is compared to the previous lot by mass spectrometry to ensure maximum lot-to-lot

consistency.



## **DATASHEET**

Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK Telephone: +44 (0) 1223 755950 - Fax: +44 (0) 1223 755951 - E-Mail: info@abbexa.com

**Note:** This product is for research use only.