

### **Biotinylated Anti-Rat Leptin Antibody**

Catalog # ABG10344

### **Specification**

## Biotinylated Anti-Rat Leptin Antibody - Product Information

Application WB, E
Reactivity Rat
Host Goat
Clonality Polyclonal

Biotinylated Anti-Rat Leptin Antibody - Additional Information

#### **Preparation**

Produced from sera of goats pre-immunized with highly pure (>98%) recombinant Rat Leptin. Anti-Rat Leptin specific antibody was purified by affinity chromatography and then biotinylated.

#### WesternBlot


To detect Rat Leptin by Western Blot analysis this antibody can be used at a concentration of 0.1 -  $0.2~\mu g/ml$ . Used in conjunction with compatible secondary reagents the detection limit for recombinant Rat Leptin is 1.5 - 3.0~ng/lane, under either reducing or non-reducing conditions.

#### Sandwich

To detect Rat Leptin by sandwich ELISA (using  $100 \, \mu \text{I/well}$  antibody solution) a concentration of  $0.25 - 1.0 \, \mu \text{g/ml}$  of this antibody is required. This biotinylated polyclonal antibody, in conjunction with BioGems' Polyclonal Anti-Rat Leptin (62-063P) as a capture antibody, allows the detection of at least  $0.2 - 0.4 \, \text{ng/well}$  of recombinant Rat Leptin.

### Direct

To detect Rat Leptin by direct ELISA (using 100  $\mu$ l/well antibody solution) a concentration of 0.25 – 1.0  $\mu$ g/ml of this antibody is required. This biotinylated polyclonal antibody, in conjunction with compatible secondary reagents, allows the detection of at least 0.2 – 0.4 ng/well of recombinant Rat Leptin.





### **Formulation**

A sterile filtered antibody solution was lyophilized from PBS, pH 7.2.

## Reconstitution

Centrifuge vial prior to opening. Reconstitute in sterile PBS containing 0.1% BSA to a concentration of 0.1-1.0 mg/ml.

## Storage -20°C

#### **Precautions**

Biotinylated Anti-Rat Leptin Antibody is for research use only and not for use in diagnostic or therapeutic procedures.

# **Biotinylated Anti-Rat Leptin Antibody - Protocols**

Provided below are standard protocols that you may find useful for product applications.

- Western Blot
- Blocking Peptides
- Dot Blot
- <u>Immunohistochemistry</u>
- <u>Immunofluorescence</u>
- Immunoprecipitation
- Flow Cytomety
- Cell Culture