


Rat Cardiac Ventricular Myocytes (Adult)
ORDER INFORMATION

Name of Cells: Rat Cardiac Ventricular Myocytes (RCVMyo-ad)
Catalogue Number: **cAP-r0009-ad**
Product Format: Frozen Vial
Cell Number: >1 x 10⁶cells/vial

General Information

RCVMyo-ad cells are isolated from normal 4-week old SD rats (From Charles River, MA USA) ventricular tissue samples and shipped in a frozen vial (> 1x10⁶cells/vial, @ passage 2). Cardiac Myocytes Growth Medium (cAP-41, containing 10% Fetal Bovine Serum, and supplement factors) is recommended for cell culture. **RCVMyo-ad** cells have limited capacity for proliferation, and we do not recommend the long term culture of **RCVMyo-ad**.



cAP-r0009ad Rat Cardiac Ventricular Myocytes (adult)
expressing alpha-actinin

Characterization: RCVMyo-ad cells are tested positive for cardiomyocyte markers including alpha smooth muscle actin, alpha actinin, desmin, and Troponin T Cardiac. RCVMyo-ad cells are tested negative for common experimental animal pathogens (screen by Charles River, MA USA) and mycoplasma in vitro.

Product Use: RCVMyo-ad cells are for Research Use Only.

Shipping: Frozen vial in Dry Ice Package.

Handling of Arriving Cells

When you receive the cells in a frozen vial, you can transfer the vial of cells into a -80°C freezer for short period storage or a liquid nitrogen tank for long term storage. Thaw the cells in a 37°C water bath, and then transfer the cells in a T75 flask pre-coated with Quick coating solution (cAP-01), containing 15ml full medium.

RCVMyo cells have limited capacity for proliferation.

1. Subculture Protocol:

- A) Rinse the cells in T75 flask with 10ml of HBSS w/o Ca²⁺ and Mg²⁺ (cAP-11) at (Room Temperature, **RT**) twice.
- B) Add 2ml of Trypsin/EDTA Solution (**RT**) (cAP-23) into T75 flask (make sure the whole surface of the T75 flask is covered with Trypsin/EDTA), and gently dispose the Trypsin/EDTA solution **within 10 seconds** with aspiration.
- C) Leave the T75 flask with the cells at **RT** for 2 minutes (the HNDFCs usually will be detached from the surface within 2 minute).
- D) Add 5ml Trypsin Neutralization Buffer and spin the cells down with 800g for 5 minutes.
- E) Re-suspend the cells with 45ml of fresh Full medium and the cell suspension is transferred directly into 3 x T75 flasks (15ml each, and the cells are sub-cultured at 1:3 ratio).
- F) Culture medium (full medium) is changed every 2-3 days.

Related Products:

Quick Coating Solution	cAP-01	240ml	Angio-Proteomie
HBSS w/o Ca ²⁺ , Mg ²⁺	cAP-11	100ml	Angio-Proteomie
Cardiac Myocytes Growth Medium	cAP-41	500ml	Angio-Proteomie
Cell Freezing Solution (FBS)	cAP-22	50ml	Angio-Proteomie
Cell Freezing Solution (Non-FBS)	cAP-22B	50ml	Angio-Proteomie
Trypsin/EDTA Solution	cAP-23	100ml	Angio-Proteomie
Trypsin Neutralization Solution	cAP-28	100ml	Angio-Proteomie
ITS (100x)	cAP-26	10ml	Angio-Proteomie
L-Glutamine-MAXIMUM (100x)	cAP-27	100ml	Angio-Proteomie
Human Plasma Fibronectin Solution	cAP-42	1mg/ml	Angio-Proteomie

THESE PRODUCTS ARE FOR RESEARCH USE ONLY

Caution: Handling human and animal tissue derived products is potentially bio-hazardous. Although each cell strain is tested negative for HIV, HBV and HCV DNA, or pathogens, diagnostic tests are not necessarily 100% accurate; therefore proper precautions must be taken to avoid inadvertent exposure. Always wear gloves and safety glasses when working with these materials. Never mouth pipette. We recommend following the universal procedures for handling products of human origin as the minimum precaution against contamination.

