

Rabbit Anti-Rat Helospectin pAb

Catalog No.	PS238	Quantity:	50 µl
Specificity:	Positive control: Stefanini-fixed frozen sections of rat small intestine. (1% carbodiimide added to the fixative). Absorption with 10-100 µg helospectin per ml diluted antiserum abolishes the staining, while secretin, PHI and PACAP do not. Cross-reacts with helodermin.		
Host:	Rabbit		
Immunogen:	Synthetic helospectin I. Tested on rat and mouse tissues.		
Formulation:	50 µl undiluted, freeze-dried antiserum. Dissolve the antiserum in 50-100 µL distilled water, and dilute further in 0.1M PBS with 1% BSA and 0.1% NaN ₃ .		
Applications:	Works on: Frozen sections. 1:400-1:800 using immunofluorescence, with over night incubation at +4°C.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

