

## CDK338 Accessory Pack for 5-Plate ELISA Set

## **ELISA Set Buffers & Materials for Optimal Performance**

MATERIAL	QUANTITY	DESCRIPTION
96-well microtiter plate	5	12 removable strips
Plastic plate covers	10	Adhesive sealers
Coated plate storage bag	5	Plastic bag and desiccant
Reconstitution Buffer (0.1 ml)	1	10X concentrate, add 0.9 ml water to vial
Coating Buffer (0.5 ml)	1	10X concentrate, sterile
Wash Buffer (25 ml)	1	200X concentrate
Blocking Buffer (10 ml)	5	Ready to use. Sterile
Standard Dilution Buffer (12.5 ml)	1	10X concentrate, for Standards, samples
Secondary Antibody Dilution Buffer (25 ml)	1	Ready to use
HRP Dilution Buffer (25 ml)	2	Ready to use.
Stop Reagent (10 ml)	5	Ready to use, 1N Sulfuric Acid

## **Storage Instructions:**

Store unopened accessory pack at 2-8°C.

All reagents should be warmed to room temperature before use. If crystals have formed in concentrated Buffer (Wash and Standard Dilution): warm it gently until they have completely dissolved.

Immediately after use remaining reagents and should be returned to cold storage (2-8°C).

Once prepared, store diluted buffers 1X at 2-8°C for up to one week.

Expiry of the reagents is stated on box front labels. Do not use components past the expiration date.

## Safety & precautions for use:

The Stop Solution provided is an acid solution.

Reconstitution, Standard Dilution, Secondary Antibody Dilution and HRP Dilution Buffers contain ProClin with may cause an allergic skin reaction.

Wear appropriate clothing, gloves and eye protection.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

Toll Free: 888 769-1246

978-572-1070

978-992-0298

Phone:

Fax:

E-mail: info@cellsciences.com

Web Site: www.cellsciences.com

