

Human IL-10 ELISPOT Kit with Sterile PVDF Plates

Catalog No: CKH010A
CKH010B

Lot: xxxxxxxx


Size: 2-plate format
5-plate format

Expiration Date: xx.xx.xxxx


Contents of the kit:

| Catalog # | Item | Storage Temp. | Quantity (2-plate format) | Quantity (5-plate format) |
|-----------|--------------------------------------------------|------------------|------------------------------|------------------------------|
| CKH010-A | Coating antibody (lyophilized) | 4°C | 1 vial | 1 vial |
| CKH010-B | Biotinylated detection antibody (lyophilized) | 4°C | 1 vial | 1 vial |
| CKH010-C | Streptavidin-HRP conjugate (lyophilized) | ≤-20°C | 1 vial | 1 vial |
| | AEC coloring system: | | | |
| CKH010-D | I. AEC stock solution | ≤-20°C (in dark) | 1 vial (4 ml) | 1 vial (4 ml) |
| CKH010-E | II. Substrate buffer capsule | RT | 1 vial (2.5 ml) | 1 vial (5 ml) |
| CKH010-F | Blocking stock solution (10x) | 4°C | 1 vial (4 ml) | 1 vial (10 ml) |
| CKH010-G | Dilution buffer R (10x) | 4°C | 1 vial (4 ml) | 1 vial (10 ml) |
| CKH010-H | Tween-20 | RT | 1 vial (5 ml) | 1 vial (5 ml) |
| CKH010-P | Sterile 96-well PVDF plate (Millipore MSIPS4510) | RT | 2 | 5 |
| CKH010-S | Adhesive cover slip | RT | 5 | 10 |

Typical data for CKH010:



Stimulus: none (4x10⁵ cells/well)



Stimulus: PMA/ionomycin (5x10⁴ cells/well)

Example of IL-10 specific spots produced by human PBMC.
Preincubation: 42h, incubation ELISPOT: 20h.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences, Inc.
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Web Site: www.cellsciences.com

cellsciences.com



Cell Sciences, Inc.
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Web Site: www.cellsciences.com