

## Mouse IL-10 ELISPOT Kit with Sterile PVDF Plates

**Catalog No:** CKM031A  
CKM031B

**Lot:** xxxxxxxx

**Size:** 2-plate format  
5-plate format

**Expiration Date:** xx.xx.xxxx


### Contents of the kit:

Catalog #	Item	Storage Temp.	Quantity (2-plate format)	Quantity (5-plate format)
CKM031-A	Coating antibody (lyophilized)	4°C	1 vial	1 vial
CKM031-B	Biotinylated detection antibody (lyophilized)	4°C	1 vial	1 vial
CKM031-C	Streptavidin-HRP conjugate (lyophilized)	≤-20°C	1 vial	1 vial
<b>AEC coloring system:</b>				
CKM031-D	I. AEC stock solution	≤-20°C (in dark)	1 vial (4 ml)	1 vial (4 ml)
CKM031-E	II. Substrate buffer capsule	RT	1 vial (2.5 ml)	1 vial (5 ml)
CKM031-F	Blocking stock solution (10x)	4°C	1 vial (4 ml)	1 vial (10 ml)
CKM031-G	Dilution buffer R (10x)	4°C	1 vial (4 ml)	1 vial (10 ml)
CKM031-H	Tween-20	RT	1 vial (5 ml)	1 vial (5 ml)
CKM031-P	Sterile 96-well PVDF plate (Millipore MSIPS4510)	RT	2	5
CKM031-S	Adhesive cover slip	RT	5	10

### Typical data for CKM031



Stimulus: none ( $2 \times 10^5$  cells/well)



Stimulus: Concanavalin A ( $1 \times 10^4$  cells/well)

Example of IL-10 specific spots produced by mouse (Balb/c) splenocytes.  
Preincubation: 24h, incubation ELISPOT: 24h.

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**



**Cell Sciences, Inc.**  
65 Parker Street  
Unit 11  
Newburyport, MA 01950

Toll Free: 888-769-1246  
Phone: 978-572-1070  
Fax: 978-992-0298

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Web Site: [www.cellsciences.com](http://www.cellsciences.com)

# cellsciences.com



**Cell Sciences, Inc.**  
65 Parker Street  
Unit 11  
Newburyport, MA 01950

Toll Free: 888-769-1246  
Phone: 978-572-1070  
Fax: 978-992-0298

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Web Site: [www.cellsciences.com](http://www.cellsciences.com)