



Mite Extract-Dp

| | |
|----------------|---|
| Form | Lyophilized protein powder Mite-Dp Crude Extract |
| | This product extracted from House Dust Mites, Dermatophagoides pteronyssinus (Trouessari)*. *Mites were homogenized with 10-fold 50 mM Phosphate buffer (pH7.2) and the homogenate were stirred at 4°C for 24 hours. And centrifugal supernatant dialyzed 4 times by PBS(Visking tube, 4°C). |
| Volume | Protein: 10mg |
| Storage | Store below 4°C (below -20°C for prolonged storage or after dissolved) |
| Note | LSL allergen extract is lyophilized product containing totally 10mg proteins per vial. Therefore, volume of distilled water to reconstitute is different by production lots. For detailed information, visit the website below https://www.cosmobio.com/products/ls_lg5449_1.html |

Allergic antigen

| | | | |
|---------|-----------------------------|---------|--------------------------------------|
| LG-0009 | DNP-Ascaris | LG-1109 | TNP-Ascaris |
| LG-0533 | DNP-Mite-Df | LG-1117 | TNP-BSA |
| LG-0544 | DNP-Mite-Dp | LG-2334 | Mite-Df Feces AG |
| LG-0528 | DNP-Cedar Pollen-Cj | LG-2444 | Mite-Dp Feces AG |
| LG-0017 | DNP-BSA | LG-5009 | Ascaris Extract |
| LG-3017 | DNP-BSA (30) | LG-5339 | Mite Extract-Df |
| LG-0024 | DNP-OA (Ovalbumin) | LG-5449 | Mite Extract-Dp |
| LG-0047 | DNP-Casein | LG-5559 | Mite Extract-Tp |
| LG-0067 | DNP-LG (BETA-Lactoglobulin) | LG-5229 | Cedar Pollen Extract |
| LG-0089 | DNP-KLH | LG-5280 | Cedar Pollen Extract-Cj |
| LG-0317 | BPO-BSA | LG-5779 | Hinoki Cypress, Pollen Crude Extract |
| LG-0577 | BPO-BGG | LG-5889 | Ragweed Pollen Extract |

Adjuvant

| | |
|---------|------|
| LG-6000 | ALUM |
|---------|------|

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp> e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617 [国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618 FAX : +81-3-5632-9619

Mite Extract-Dp

I. 内容

lot No. 113031

| | |
|-----|---|
| 性 状 | 凍結乾燥粉末 (Mite-Dp粗抽出物) House dust mites, <i>Dermatophagoides pteronyssinus</i> (Trouessari) ヤケヒョウヒダニ (チリダニ科) の成虫から抽出* |
| 規 格 | 蛋白質質量: 10mg 3.6mlの精製水を加える事によりPBS (Dulbecco, pH7.4) 溶液になります。 |


* 虫体を10倍量の50mMリン酸緩衝液 (pH7.2) でホモジェナイズ、24時間攪拌 (4℃)、遠心上澄を濃縮後PBSで透析 (Visking tube, 4℃)。 既知アレルゲンの含有量はGroup1が約80μg/vial、Group2が約300μg/vial程度です。

II. 保存上の注意

4℃以下 (長期間又は溶液の場合は-20℃以下) で保存して下さい。

製造元

総発売元

 株式会社 エル・エス・エル


コスモ・バイオ株式会社

LSL社の抽出物は、蛋白質抽出時の濃度がロットによって異なります。ロット毎に加える精製水の量も異なりますので、商品に添付されているデータシートをご確認ください。

新製品案内

| | |
|---------|---------------------------------|
| LG-2334 | Mite-Df Feces AG (ダニーDf由来排泄物抗原) |
| LG-2444 | Mite-Dp Feces AG (ダニーDp由来排泄物抗原) |
| LG-5999 | Wormwood Pollen Extract |
| LG-7000 | FTPP Suspension (百日咳死菌懸濁液) |
| LG-8000 | FTMB Suspension (結核死菌懸濁液) |
| LB-7111 | Anti (MAP) Mite Der fl |
| LB-7103 | Anti (MAP) Mite Group 1 |
| LB-7204 | Anti (MAP) Mite Group 2 |