

Actin Polymerization Buffer

V. 2.0

Cat. # BSA02-001

Lot # Amount : 1 x 2 ml

Upon arrival store at 4°C (desiccated)

See datasheet for storage after reconstitution

Material

Actin Polymerization Buffer is supplied as a white lyophilized powder.

Storage and Reconstitution

Lyophilized Actin Polymerization Buffer is stable at 4°C desiccated (<10% humidity) for 1 year.

Each tube of Actin Polymerization Buffer should be resuspended in 2 ml of 100 mM Tris HCl pH 7.5 to give a 10X strength buffer; 100 mM Tris HCl, 20 mM MgCl₂, 500 mM KCl, 10 mM ATP, 50 mM guanidine carbonate pH 7.5. The reconstituted buffer should be aliquoted into experiment sized volumes, snap frozen in liquid nitrogen and stored at -70°C.

Uses

A frozen aliquot of the Actin Polymerization Buffer should be thawed at room temperature and placed on ice immediately after thawing. The buffer contains ATP which is a labile component. Hence, the buffer should be used within 1-2 h of thawing and freeze/thaw cycles should be avoided.

Tris HCl is present in the buffer to maintain a pH of 7.5. Magnesium chloride (2 mM final concentration), potassium chloride (50 mM final concentration) and ATP (1 mM final concentration) are required for actin polymerization. Guanidine carbonate is present as a stabilizer for ATP and is not required for actin polymerization.

Used at 1X final concentration as an *in vitro* buffer for monomer actin polymerization from a wide range of species.

Product Citations / Related Products

For the latest citations and related products please visit www.cytoskeleton.com.