Rhesus CD3d/CD3 delta Gene ORF cDNA clone expression plasmid, C-GFPSpark tag



Catalog Number: CG90046-ACG

General Information

Gene :	CD3d molecule, delta (CD3-TCR complex)	
Official Symbol :	CD3D	
Synonym :	CD3D	
Source :	Rhesus	
cDNA Size:	516bp (cDNA Size= Gene + linker +Tags)	
RefSeq :	XM_015115818.2	
Plasmid:	pCMV3-cynoCD3D-GFPSpark	
-		

Description

Lot : Please refer to the label on the tube

Sequence Description :

Identical with the Gene Bank Ref. ID sequence. Please check the sequence information before order.

Restriction site: Kpnl + Xbal (6kb + 1.25kb)

Vector : pCMV3-C-GFPSpark

Quality control :

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

Sequencing primer list :

pCMV3-F:	5' CAGGTGTCCACTCCCAGGTCCAAG 3'
pcDNA3-R :	5' GGCAACTAGAAGGCACAGTCGAGG 3'
Or	
Forward T7 :	5' TAATACGACTCACTATAGGG 3'
ReverseBGH :	5' TAGAAGGCACAGTCGAGG 3'

pCMV3-F and pcDNA3-R are designed by Sino Biological Inc. Customers can order the primer pair from any oligonucleotide supplier.

Shipping carrier :

Each tube contains approximately 10 µg of lyophilized plasmid. **Storage :**

The lyophilized plasmid can be stored at ambient temperature for three months.

Plasmid Resuspension protocol

- 1. Centrifuge at 5,000 \times g for 5 min.
- 2. Carefully open the tube and add 100 μl of sterile water to

dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room

temperature.

4. Briefly vortex the tube and then do a quick spin to concentrate

the liquid at the bottom. Speed is less than 5000 imes g.

5. Store the plasmid at $\$ -20 $^\circ\! {\rm C}.$

The plasmid is ready for:

- Restriction enzyme digestion
- PCR amplification
- E. coli transformation
- DNA sequencing

E.coli strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5 α and TOP10F['].</sup>

Rhesus CD3d/CD3 delta Gene ORF cDNA clone expression plasmid, C-GFPSpark tag



Catalog Number: CG90046-ACG

Vector Information

All of the pCMV vectors are designed for high-level stable and transient expression in mammalian hosts. High-level stable and non-replicative transient expression can be carried out in most mammalian cells. The vectors contain the following elements:

•Human enhanced cytomegalovirus immediate-early (CMV) promoter for high-level expression in a wide range of mammalian cells.

• Hygromycin resistance gene for selection of mammalian cell lines.

• A Kozak consensus sequence to enhance mammalian expression.

Vector name	pCMV3-C-GFPSpark
Vector size	6848bp
Vector Type	Mammalian Expression Vector
Expression Method	Constiutive ,Stable / Transient
Promoter	CMV
Bacterial Resistance	Kanamycin
Selection In Cells	Hygromycin
Protein tag	GFPSpark

Physical Map of Plasmid :

