## Human MERTK/Mer Gene ORF cDNA clone expression plasmid, C-GFPSpark tag



**Catalog Number:** HG10298-ACG

**General Information** 

Gene: c-mer proto-oncogene tyrosine kinase

Official Symbol: **MERTK** 

Synonym: c-Eyk; c-mer; MER; RP38; Tyro12

Source: Human

cDNA Size: 3729bp

RefSeq: NM\_006343.2

Plasmid: pCMV3-MERTK-GFPSpark

**Description** 

Please refer to the label on the tube Lot:

**Sequence Description:** 

Identical with the Gene Bank Ref. ID sequence.

**Restriction site:** HindIII + Xbal(6kb+3.73kb)

Vector: pCMV3-C-GFPSpark

**Quality control:** 

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

Sequencing primer list:

pCMV3-F: 5' CAGGTGTCCACTCCCAGGTCCAAG 3' 5' GGCAACTAGAAGGCACAGTCGAGG 3' pcDNA3-R: Or

Forward T7:

5' TAATACGACTCACTATAGGG 3'

ReverseBGH: 5' TAGAAGGCACAGTCGAGG 3'

pCMV3-F and pcDNA3-R are designed by Sino Biological Inc. Customers can order the primer pair from any oligonucleotide supplier.

#### Shipping carrier:

Each tube contains approximately 10 µg of lyophilized plasmid.

#### Storage:

The lyophilized plasmid can be stored at ambient temperature for three months.

#### **Plasmid Resuspension protocol**

- 1. Centrifuge at 5,000×g for 5 min.
- 2. Carefully open the tube and add 100  $\mu l$  of sterile water to dissolve the DNA.
- 3. Close the tube and incubate for 10 minutes at room temperature.
- 4. Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than  $5000 \times g$ .
- 5. Store the plasmid at -20 °C.

#### The plasmid is ready for:

- Restriction enzyme digestion
- · PCR amplification
- · E. coli transformation
- · DNA sequencing

### E.coli strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5α and TOP10F'.

Fax:+86-10-51029969

# Human MERTK/Mer Gene ORF cDNA clone expression plasmid, C-GFPSpark tag



Catalog Number: HG10298-ACG

#### **Vector Information**

All of the pCMV vectors are designed for high-level stable and transient expression in mammalian hosts. High-level stable and non-replicative transient expression can be carried out in most mammalian cells. The vectors contain the following elements:

- •Human enhanced cytomegalovirus immediate-early (CMV) promoter for high-level expression in a wide range of mammalian cells.
- Hygromycin resistance gene for selection of mammalian cell lines.
- A Kozak consensus sequence to enhance mammalian expression.

Vector Name pCMV3-C-GFPSpark

Vector Size 6848bp

Vector Type Mammalian Expression Vector

Expression Method Constitutive, Stable / Transient

Promoter CMV

Antibiotic Resistance Kanamycin

Selection In Mammalian Cells

n Cells Hygromycin

Protein Tag GFPSpark

#### **Physical Map of Plasmid:**

