Human ICAM-1/CD54 Gene ORF cDNA clone expression plasmid, C-His tag

5' TAATACGACTCACTATAGGG 3'

5' TAGAAGGCACAGTCGAGG 3'

pCMV3-F and pcDNA3-R are designed by Sino Biological Inc. Customers can order the primer pair from any oligonucleotide supplier.



Catalog Number: HG10346-CH

Or

Forward T7 :

ReverseBGH :

General Information		Shipping carrier :
Gene :	intercellular adhesion molecule 1	Each tube contains approximately 10 µg of lyophilized plasmid.
Gene .		Storage :
Official Symbol :	ICAM1	The lyophilized plasmid can be stored at ambient temperature for
Synonym :	BB2; CD54; ICAM-1; P3.58	three months.
Source :	Human	
cDNA Size:	1644bp	Plasmid Resuspension protocol
RefSeq :	NM_000201.2	1. Centrifuge at 5,000 $ imes$ g for 5 min.
Plasmid:	pCMV3-ICAM1-His	2. Carefully open the tube and add 100 μI of sterile water to
Description		dissolve the DNA.
Lot: Please	refer to the label on the tube	3. Close the tube and incubate for 10 minutes at room
Sequence Description :		temperature.
Identical with the Gene Bank Ref. ID sequence.		4. Briefly vortex the tube and then do a quick spin to concentrate
Restriction site:	Kpnl + Xbal(6kb+1.64kb)	the liquid at the bottom. Speed is less than 5000 $ imes$ g.
Vector :	pCMV3-C-His	5. Store the plasmid at -20 °C.
Quality control :		
The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.		The plasmid is ready for:
Sequencing primer list :		Restriction enzyme digestion
pCMV3-F:	5' CAGGTGTCCACTCCCAGGTCCAAG 3'	PCR amplification
pcDNA3-R :	5' GGCAACTAGAAGGCACAGTCGAGG 3'	• E. coli transformation

DNA sequencing

E.coli strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5 α and TOP10F'.

Human ICAM-1/CD54 Gene ORF cDNA clone expression plasmid, C-His tag



Catalog Number: HG10346-CH

Vector Information


All of the pCMV vectors are designed for high-level stable and transient expression in mammalian hosts. High-level stable and non-replicative transient expression can be carried out in most mammalian cells. The vectors contain the following elements:

•Human enhanced cytomegalovirus immediate-early (CMV) promoter for high-level expression in a wide range of mammalian cells.

• Hygromycin resistance gene for selection of mammalian cell lines.

• A Kozak consensus sequence to enhance mammalian expression.

Physical Map of Plasmid :



Vector Name	pCMV3-C-His
Vector Size	6164bp
Vector Type	Mammalian Expression Vector
Expression Method	Constitutive, Stable / Transient
Promoter	CMV
Antibiotic Resistance	Kanamycin
Selection In Mammalian Cells	Hygromycin
Protein Tag	His