Human CCL4 ORF mammalian expression plasmid, N-HA tag



Catalog Number: HG10899-NY

General Information

Gene: chemokine (C-C motif) ligand 4

Official Symbol: CCL4

Synonym: ACT2, G-26, LAG1, MIP1B, SCYA2,

> SCYA4, MIP1B1, AT744.1, MGC104418, MGC126025, MGC126026, MIP-1-beta

Source: Human

cDNA Size: 279bp

RefSeq: NM_002984.2

Plasmid: pCMV3-HA-CCL4

Description

Lot: Please refer to the label on the tube

Sequence Description:

Identical with the Gene Bank Ref. ID sequence except for the point mutations: 138T/C not causing the amino acid variation.

Restriction site: Kpnl + Xbal (6kb + 0.3kb)

Vector: pCMV3-SP-N-HA

Shipping carrier:

Each tube contains approximately 10 µg of lyophilized plasmid.

Storage:

The lyophilized plasmid can be stored at ambient temperature for three months.

Quality control:

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

Sequencing primer list:

pCMV3-F: 5' CAGGTGTCCACTCCCAGGTCCAAG 3'

pcDNA3-R: 5' GGCAACTAGAAGGCACAGTCGAGG 3'

Or

Forward T7: 5' TAATACGACTCACTATAGGG 3'

ReverseBGH: 5' TAGAAGGCACAGTCGAGG 3'

pCMV3-F and pcDNA3-R are designed by Sino Biological Inc. Customers can order the primer pair from any oligonucleotide supplier.

Plasmid Resuspension protocol

- 1. Centrifuge at 5,000×g for 5 min.
- Carefully open the tube and add 100 µl of sterile water to dissolve the DNA.
- 3. Close the tube and incubate for 10 minutes at room temperature.
- Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than 5000×g.
- 5. Store the plasmid at -20 °C.

The plasmid is ready for:

- · Restriction enzyme digestion
- PCR amplification
- · E. coli transformation
- DNA sequencing

E.coli strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5α and TOP10F'.

Human CCL4 ORF mammalian expression plasmid, N-HA tag



Catalog Number: HG10899-NY

Vector Information

All of the pCMV vectors are designed for high-level stable and transient expression in mammalian hosts. High-level stable and non-replicative transient expression can be carried out in most mammalian cells. The vectors contain the following elements:

- •Human enhanced cytomegalovirus immediate-early (CMV) promoter for high-level expression in a wide range of mammalian cells.
- Hygromycin resistance gene for selection of mammalian cell lines.
- A Kozak consensus sequence to enhance mammalian expression.

Vector Name pCMV3-SP-N-HA

Vector Size 6146bp

Vector Type Mammalian Expression Vector

Expression Method Constitutive, Stable / Transient


Promoter CMV

Antibiotic Resistance Kanamycin

Selection In Mammalian Cells Hygromycin


Protein Tag HA

Physical Map of Plasmid:



pCMV3-SP-N-HA (suitable for secretary and membane protein expession)





Comments for pCMV3-SP-N-HA:


CMV promoter: bases 250-837 enhancer: bases 838-1445

SV40 early promoter: bases 2387-2756 Hygromycin ORF: bases 2774-3799 pUC origin: bases 4442-5115 Kanamycin ORF: bases 5189-6004

Description

Vector Name	pCMV3-SP-N-HA
Vector Size	6146bp
Vector Type	Mammalian Expression Vector
Expression Method	Constitutive, Stable / Transient
Promoter	CMV
Antibiotic Resistance	Kanamycin
Selection In Mammalian Cells	Hygromycin
Protein Tag	HA
Sequencing Primer	Forward:T7(TAATACGACTCACTATAGGG) Reverse:BGH(TAGAAGGCACAGTCGAGG)

Schematic of pCMV3-SP-N-HA Multiple Cloning Sites



pCMV3-SP-N-HA is recommended for constructing the N-HA tag secretary and membrane proteins expression vector which containing a naïve signal peptide. An universal signal peptide is used to instead the naïve signal peptide.