# Human SAA1 Gene ORF cDNA clone expression plasmid, C-GFPSpark tag



Catalog Number: HG11298-ACG

**General Information** 

**Gene:** serum amyloid A1

Official Symbol: SAA1

Synonym: PIG4; SAA; SAA2; TP53I4

Source: Human

cDNA Size: 1098bp

**RefSeq:** NM\_000331.3

Plasmid: pCMV3-SAA1-GFPSpark

Description

Lot: Please refer to the label on the tube

**Sequence Description:** 

Identical with the Gene Bank Ref. ID sequence.

**Restriction site:** KpnI + XbaI(6kb+1.1kb)

**Vector**: pCMV3-C-GFPSpark

**Quality control:** 

Forward T7:

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

Sequencing primer list:

pCMV3-F: 5' CAGGTGTCCACTCCCAGGTCCAAG 3'
pcDNA3-R: 5' GGCAACTAGAAGGCACAGTCGAGG 3'
Or

5' TAATACGACTCACTATAGGG 3'

ReverseBGH: 5' TAGAAGGCACAGTCGAGG 3'

pCMV3-F and pcDNA3-R are designed by Sino Biological Inc. Customers can order the primer pair from any oligonucleotide supplier.

#### Shipping carrier:

Each tube contains approximately 10 µg of lyophilized plasmid.

#### Storage:

The lyophilized plasmid can be stored at ambient temperature for three months.

### **Plasmid Resuspension protocol**

- 1. Centrifuge at 5,000×g for 5 min.
- 2. Carefully open the tube and add 100  $\mu l$  of sterile water to dissolve the DNA.
- 3. Close the tube and incubate for 10 minutes at room temperature.
- 4. Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than  $5000 \times g$ .
- 5. Store the plasmid at  $\ -20\ ^{\circ}\mathrm{C}$  .

#### The plasmid is ready for:

- Restriction enzyme digestion
- PCR amplification
- E. coli transformation
- · DNA sequencing

### *E.coli* strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5 $\alpha$  and TOP10F $^{\prime}$ .

# Human SAA1 Gene ORF cDNA clone expression plasmid, C-GFPSpark tag



Catalog Number: HG11298-ACG

#### **Vector Information**

All of the pCMV vectors are designed for high-level stable and transient expression in mammalian hosts. High-level stable and non-replicative transient expression can be carried out in most mammalian cells. The vectors contain the following elements:

- •Human enhanced cytomegalovirus immediate-early (CMV) promoter for high-level expression in a wide range of mammalian cells.
- Hygromycin resistance gene for selection of mammalian cell lines.
- A Kozak consensus sequence to enhance mammalian expression.

Vector Name pCMV3-C-GFPSpark

Vector Size 6848bp

Vector Type Mammalian Expression Vector

Expression Method Constitutive, Stable / Transient

Promoter CMV

Antibiotic Resistance Kanamycin

Selection In Mammalian Cells

Hygromycin

Protein Tag

**GFPSpark** 

#### **Physical Map of Plasmid:**

