Human MAPK3 ORF mammalian expression plasmid, N-OFPSpark / RFP tag



Catalog Number: HG11545-ANR

General Information

Gene : mitogen-activated protein kinase 3

Official Symbol :	MAPK3
Synonym :	ERK1, PRKM3, P44ERK1, P44MAPK, HS44KDAP, HUMKER1A, MGC20180, MAPK3
Source :	Human

RefSeq : NM_002746.2

Description

cDNA Size:

Lot : Please refer to the label on the tube

1140bp

Vector : pCMV3-N-OFPSpark

Shipping carrier :

Each tube contains approximately 10 µg of lyophilized plasmid.

Storage :

The lyophilized plasmid can be stored at ambient temperature for three months.

Quality control :

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

Sequencing primer list :

pCMV3-F:	5' CAGGTGTCCACTCCCAGGTCCAAG 3'	
pcDNA3-R :	5' GGCAACTAGAAGGCACAGTCGAGG 3'	
Or		
Forward T7 :	5' TAATACGACTCACTATAGGG 3'	
ReverseBGH :	5' TAGAAGGCACAGTCGAGG 3'	
nCMV/2 E and no DNA 2 B and designed by Gine Dislogical las		

pCMV3-F and pcDNA3-R are designed by Sino Biological Inc. Customers can order the primer pair from any oligonucleotide supplier.

Plasmid Resuspension protocol

- 1. Centrifuge at $5,000 \times g$ for 5 min.
- 2. Carefully open the tube and add 100 μ l of sterile water to dissolve the DNA.
- Close the tube and incubate for 10 minutes at room temperature.
- 4. Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than $5000 \times g$.
- 5. Store the plasmid at -20 $^\circ\!\mathrm{C}.$

The plasmid is ready for:

- Restriction enzyme digestion
- PCR amplification
- E. coli transformation
- DNA sequencing

E.coli strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5 α and TOP10F $\dot{}$.

Human MAPK3 ORF mammalian expression plasmid, N-OFPSpark / RFP tag



Catalog Number: HG11545-ANR

Vector Information

All of the pCMV vectors are designed for high-level stable and transient expression in mammalian hosts. High-level stable and non-replicative transient expression can be carried out in most mammalian cells. The vectors contain the following elements:

•Human enhanced cytomegalovirus immediate-early (CMV) promoter for high-level expression in a wide range of mammalian cells.

• Hygromycin resistance gene for selection of mammalian cell lines.

• A Kozak consensus sequence to enhance mammalian expression.

Vector Name	pCMV3-N-OFPSpark
Vector Size	6746bp
Vector Type	Mammalian Expression Vector
Expression Method	Constiutive ,Stable / Transient
Promoter	CMV
Antibiotic Resistance	Kanamycin
Selection In Mammalian Cells	Hygromycin
Protein Tag	OFPSpark
Sequencing Primer	Forward:T7(TAATACGACTCACTATAGGG) Reverse:BGH(TAGAAGGCACAGTCGAGG)