## Human AR ORF mammalian expression plasmid, C-Flag tag



Catalog Number: HG13116-CF

**General Information** 

Gene: androgen receptor

Official Symbol : AR

**Synonym:** KD, AIS, TFM, DHTR, SBMA, HYSP1,

NR3C4, SMAX1, HUMARA, AR

Source: Human

cDNA Size: 1167bp

**RefSeq:** NM\_001011645.1

Plasmid: pCMV3-AR-Flag

Description

Lot: Please refer to the label on the tube

**Sequence Description:** 

Identical with the Gene Bank Ref. ID sequence.

Restriction site: Kpnl + Xbal (6kb + 1.21kb)

**Vector:** pCMV3-C-FLAG

Shipping carrier:

Each tube contains approximately 10 µg of lyophilized plasmid.

Storage:

The lyophilized plasmid can be stored at ambient temperature for three months.

Quality control:

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

Sequencing primer list:

pCMV3-F: 5' CAGGTGTCCACTCCCAGGTCCAAG 3'

pcDNA3-R: 5' GGCAACTAGAAGGCACAGTCGAGG 3'

Or

Forward T7: 5' TAATACGACTCACTATAGGG 3'

ReverseBGH: 5' TAGAAGGCACAGTCGAGG 3'

pCMV3-F and pcDNA3-R are designed by Sino Biological Inc. Customers can order the primer pair from any oligonucleotide supplier.

### **Plasmid Resuspension protocol**

- 1. Centrifuge at  $5,000 \times g$  for 5 min.
- 2.Carefully open the tube and add 100  $\mu l$  of sterile water to dissolve the DNA.
- 3. Close the tube and incubate for 10 minutes at room temperature.
- 4.Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than  $5000\times g$ .
- 5.Store the plasmid at -20 °C.

#### The plasmid is ready for:

- Restriction enzyme digestion
- · PCR amplification
- E. coli transformation
- DNA sequencing

### *E.coli* strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5 $\alpha$  and TOP10F'.

# **Human AR ORF mammalian expression** plasmid, C-Flag tag



Catalog Number: HG13116-CF

### **Vector Information**

All of the pCMV vectors are designed for high-level stable and transient expression in mammalian hosts. High-level stable and non-replicative transient expression can be carried out in most mammalian cells. The vectors contain the following elements:

- •Human enhanced cytomegalovirus immediate-early (CMV) promoter for high-level expression in a wide range of mammalian cells.
- Hygromycin resistance gene for selection of mammalian cell lines.
- A Kozak consensus sequence to enhance mammalian expression.

Vector Name pCMV3-C-FLAG

Vector Size 6158bp

Vector Type Mammalian Expression Vector Expression Method Constitutive, Stable / Transient

Promoter CMV Antibiotic Resistance Kana

Antibiotic Resistance Kanamycin Selection In Mammalian Cells Hygromycin

Protein Tag FLAG

### Physical Map of Plasmid:

