Human HSPA9/Grp75 Gene ORF cDNA clone expression plasmid, N-Flag tag



Catalog Number: HG16926-NF

General Information

heat shock 70kDa protein 9 (mortalin)
HSPA9
CRP40; CSA; GRP-75; GRP75; HEL-S- 124m; HSPA9B; MOT; MOT2; MTHSP75; PBP74
Human
2079bp (cDNA Size= Gene + linker +Tags)
NM_004134.6
pCMV3-Flag-HSPA9

Description

Lot : Please refer to the label on the tube

Sequence Description :

Identical with the Gene Bank Ref. ID sequence except for the point mutations: 1972C/T not causing the amino acid variation.

Restriction site: Kpnl + Notl(6kb+2.08kb)

Vector : pCMV3-N-FLAG

Quality control :

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

Sequencing primer list :

pCMV3-F:	5' CAGGTGTCCACTCCCAGGTCCAAG 3'
pcDNA3-R :	5' GGCAACTAGAAGGCACAGTCGAGG 3'
Or	
Forward T7 :	5' TAATACGACTCACTATAGGG 3'
ReverseBGH :	5' TAGAAGGCACAGTCGAGG 3'

pCMV3-F and pcDNA3-R are designed by Sino Biological Inc. Customers can order the primer pair from any oligonucleotide supplier.

Shipping carrier :

Each tube contains approximately 10 µg of lyophilized plasmid. **Storage :**

The lyophilized plasmid can be stored at ambient temperature for three months.

Plasmid Resuspension protocol

- 1. Centrifuge at 5,000 \times g for 5 min.
- 2. Carefully open the tube and add 100 μl of sterile water to

dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room

temperature.

4. Briefly vortex the tube and then do a quick spin to concentrate

the liquid at the bottom. Speed is less than $5000 \times g$.

5. Store the plasmid at -20 °C.

The plasmid is ready for:

- Restriction enzyme digestion
- PCR amplification
- E. coli transformation
- DNA sequencing

E.coli strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5 α and TOP10F['].</sup>

Human HSPA9/Grp75 Gene ORF cDNA clone expression plasmid, N-Flag tag



Catalog Number: HG16926-NF

Vector Information

All of the pCMV vectors are designed for high-level stable and transient expression in mammalian hosts. High-level stable and non-replicative transient expression can be carried out in most mammalian cells. The vectors contain the following elements:

•Human enhanced cytomegalovirus immediate-early (CMV) promoter for high-level expression in a wide range of mammalian cells.

• Hygromycin resistance gene for selection of mammalian cell lines.

• A Kozak consensus sequence to enhance mammalian expression.

Vector name	pCMV3-N-FLAG
Vector size	6098bp
Vector Type	Mammalian Expression Vector
Expression Method	Constiutive ,Stable / Transient
Promoter	CMV
Bacterial Resistance	Kanamycin
Selection In Cells	Hygromycin
Protein tag	FLAG

Physical Map of Plasmid :

