# Human TPX2 ORF mammalian expression plasmid, C-OFPSpark / RFP tag



#### Catalog Number: HG18711-ACR

#### **General Information**

Gene : TPX2, microtubule-associated

Official	Symbol :	TPX2	
Univia	Oynibol .	• • • • • • • • • • • • • • • • • • • •	

- Synonym :
   DIL2; p100; DIL-2; HCTP4; FLS353; HCA519; REPP86; C20orf1; C20orf2; GD:C20orf1

   Source :
   Human
- cDNA Size: 2244bp
- RefSeq : NM\_012112.4

### Description

Lot : Please refer to the label on the tube

Vector : pCMV3-C-OFPSpark

#### Shipping carrier :

Each tube contains approximately 10 µg of lyophilized plasmid.

#### Storage :

The lyophilized plasmid can be stored at ambient temperature for three months.

#### **Quality control:**

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

#### Sequencing primer list :

pCMV3-F:	5' CAGGTGTCCACTCCCAGGTCCAAG 3'	
pcDNA3-R :	5' GGCAACTAGAAGGCACAGTCGAGG 3'	
Or		
Forward T7 :	5' TAATACGACTCACTATAGGG 3'	
ReverseBGH :	5' TAGAAGGCACAGTCGAGG 3'	

pCMV3-F and pcDNA3-R are designed by Sino Biological Inc. Customers can order the primer pair from any oligonucleotide supplier.

## **Plasmid Resuspension protocol**

- 1. Centrifuge at  $5,000 \times g$  for 5 min.
- 2. Carefully open the tube and add 100  $\mu$ l of sterile water to dissolve the DNA.
- Close the tube and incubate for 10 minutes at room temperature.
- 4. Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than  $5000 \times g$ .
- 5. Store the plasmid at -20  $^{\circ}$ C.

#### The plasmid is ready for:

- Restriction enzyme digestion
- PCR amplification
- E. coli transformation
- DNA sequencing

# *E.coli* strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5 $\alpha$  and TOP10F  $\dot{}$  .

# Human TPX2 ORF mammalian expression plasmid, C-OFPSpark / RFP tag



Catalog Number: HG18711-ACR

## **Vector Information**

All of the pCMV vectors are designed for high-level stable and transient expression in mammalian hosts. High-level stable and non-replicative transient expression can be carried out in most mammalian cells. The vectors contain the following elements:

•Human enhanced cytomegalovirus immediate-early (CMV) promoter for high-level expression in a wide range of mammalian cells.

• Hygromycin resistance gene for selection of mammalian cell lines.

• A Kozak consensus sequence to enhance mammalian expression.

Vector Name	pCMV3-C-OFPSpark
Vector Size	6806bp
Vector Type	Mammalian Expression Vector
Expression Method	Constiutive ,Stable / Transient
Promoter	CMV
Antibiotic Resistance	Kanamycin
Selection In Mammalian Cells	Hygromycin
Protein Tag	OFPSpark
Sequencing Primer	Forward:T7(TAATACGACTCACTATAGGG) Reverse:BGH(TAGAAGGCACAGTCGAGG)