# Mouse ATF4 Gene cDNA clone plasmid



Catalog Number: MG52828-G

## **General Information**

Gene : activating transcription factor 4

Official Symbol : ATF4

Synonym : Atf-4, C/ATF, CREB2, TAXREB67

Source : Mouse

cDNA Size: 1050bp

**RefSeq :** NM\_009716.3

Plasmid: pGEM-mATF4

# Description

Lot : Please refer to the label on the tube

### Sequence Description :

Identical with the Gene Bank Ref. ID sequence except for the point mutation: 672A>G not causing the amino acid variation.

Vector :

pGEM-T

Shipping carrier :

Each tube contains approximately 10 µg of lyophilized plasmid.

### Storage :

The lyophilized plasmid can be stored at ambient temperature for three months.

### **Quality control :**

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

## Sequencing primer list :

M13-47 : 5' GCCAGGGTTTTCCCAGTCACGAC 3'

RV-M : 5' GAGCGGATAACAATTTCACACAGG 3'

Other M13 primers can also be used as sequencing primers.

# **Plasmid Resuspension protocol**

1.Centrifuge at  $5,000 \times g$  for 5 min.

2.Carefully open the tube and add 100  $\mu l$  of sterile water to dissolve the DNA.

3.Close the tube and incubate for 10 minutes at room temperature. 4.Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than  $5000 \times g$ . 5.Store the plasmid at -20 °C.

### The plasmid is ready for:

- Restriction enzyme digestion
- PCR amplification
- E. coli transformation
- DNA sequencing

# *E.coli* strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5 $\alpha$  and TOP10F  $\dot{}$  .


# Mouse ATF4 Gene cDNA clone plasmid

Catalog Number: MG52828-G

# **Vector Information**

The pGEM-T vector is a high-efficiency TA cloning vector which contains multiple cloning sites as shown below. The pGEM-T vector is 3.0kb in size and contains the amplicin resistance gene for selection. The coding sequence was inserted by TA cloning.

## Physical Map of pGEM-T :



For Research Use Only. Not for use in diagnostic or therapeutic procedures. Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe) Website: <u>http://www.sinobiological.com</u>

Sino Biological

Biological Solution Specialist