

# Human Contactin 5 / CNTN5 ELISA Pair Set

## Catalog Number : SEK10513

To achieve the best assay results, this manual must be read carefully before using this product and the assay is run as summarized in the General ELISA protocol.

#### BACKGROUND

Contactins are a subgroup of molecules belonging to the immunoglobulin superfamily that are expressed mainly in the nervous system. The subgroup consists of six members: Contactin-1, Contactin-2(TAG-1), Contactin-3(BIG-1), BIG-2, Contactin-5(NB-2) and NB-3. Since their identification in the late 1980s, Contactin-1 and Contactin-2 have been studied extensively. Axonal expression and the neurite extension activity of Contactin-1 and Contactin-2 attracted researchers to study the function of these molecules in axon guidance during development. Contactin-1 and Contactin-2 have come to be known as the principal molecules in the function and maintenance of myelinated neurons. In contrast, the function of the other four members of this subgroup remained unknown until recently. Contactin-5, also known as NB-2, is one of the neural recognition molecules in the contactin subgroup. Contactin-5 is expressed in brain and kidney and at very low level in placenta. In brain, Contactin-5 is highly expressed in the occipital lobe, amygdala, cerebral cortex, frontal lobe, thalamus and temporal lobe. Mice deficient in the Contactin-5 gene exhibit aberrant responses to acoustic stimuli. Contactin-5 may play a role in maturation of glutamatergic synapses in the brainstem during the final stages of auditory development. Contactin-5 gene may contribute to human neurological disorders.

### PRINCIPLE OF THE TEST


The Sino Biological ELISA Pair Set is a solid phase sandwich ELISA (Enzyme-Linked Immunosorbent Assay). It utilizes a monoclonal antibody specific for Contactin 5 / CNTN5 coated on a 96-well plate. Standards and samples are added to the wells, and any Contactin 5 / CNTN5 present binds to the immobilized antibody. The wells are washed and a horseradish peroxidase conjugated mouse anti-Contactin 5 / CNTN5 monoclonal antibody is then added, producing an antibody-antigen-antibody "sandwich". The wells are again washed and TMB substrate solution is loaded, which produces color in proportion to the amount of Contactin 5 / CNTN5 present in the sample. To end the enzyme reaction, the stop solution is added and absorbances of the microwell are read at 450 nm.

#### **INTENDED USE**

◆The Human Contactin 5 / CNTN5 ELISA Pair Set is for the quantitative determination of Human Contactin 5 / CNTN5.

•This ELISA Pair Set contains the basic components required for the development of sandwich ELISAs.

#### ASSAY PROCEDURE SUMMARY



This Pair Set has been configured for research use only and is not to be used in diagnostic procedures.

## MATERIALS PROVIDED

Bring all reagents to room temperature before use.

**Capture Antibody** – 0.3 mg/mL of rabbit anti-CNTN5 monoclonal antibody. Dilute to a working concentration of 2.0 µg/mL in CBS before coating. (Catalog: # 10513-R024)

**Detection Antibody** – 0.5 mg/mL mouse anti-CNTN5 monoclonal antibody conjugated to horseradish-peroxidase (HRP). Dilute to working concentration of 0.25  $\mu$ g/mL in detection antibody dilution buffer before use. (Catalog: # 10513-MM11)

**Standard** – Each vial contains 200 ng of recombinant CNTN5. Reconstitute standard powder with 1 mL detection antibody dilution buffer. After reconstitution, store at  $-20^{\circ}$ C to  $-80^{\circ}$ C in a manual defrost freezer. A seven-point standard curve using 2-fold serial dilutions in sample dilution buffer, and a high standard of 4000 pg/mL is recommended.

## SOLUTIONS REQUIRED

CBS - 0.05M Na<sub>2</sub>CO<sub>3</sub> , 0.05M NaHCO<sub>3</sub> , pH 9.6, 0.2 µm filtered TBS - 20 mM Tris, 150 mM NaCl, pH 7.4 Wash Buffer - 0.05% Tween20 in TBS, pH 7.2 - 7.4 Blocking Buffer - 2% BSA in Wash Buffer Sample dilution buffer - 0.1% BSA in wash buffer, pH 7.2 - 7.4, 0.2 µm filtered Detection antibody dilution buffer - 0.5% BSA in wash buffer, pH 7.2 - 7.4, 0.2 µm filtered Substrate Solution : To achieve best assay results, fresh substrate solution is recommended Substrate stock solution - 10mg / ml TMB ( Tetramethylbenzidine ) in DMSO Substrate dilution buffer - 0.05M Na<sub>2</sub>HPO<sub>4</sub> and 0.025M citric acid ; adjust pH to 5.5 Substrate working solution - For each plate dilute 250 µl substrate stock solution in 25ml substrate dilution buffer and then add 80 µl 0.75% H<sub>2</sub>O<sub>2</sub> , mix it well Stop Solution - 2 N H<sub>2</sub>SO<sub>4</sub>

## PRECAUTION

The Stop Solution suggested for use with this Pair Set is an acid solution. Wear eye, hand, face, and clothing protection when using this material.

## STORAGE

**Capture Antibody**: Aliquot and store at -20 $^{\circ}$ C to -80 $^{\circ}$ C for up to 6 months from date of receipt. Avoid repeated freeze-thaw cycles.

**Detection Antibody**: Protect it from prolonged exposure to light. Aliquot and store at  $-20^{\circ}$ C to  $-80^{\circ}$ C and for up to 6 months from date of receipt. Avoid repeated freeze-thaw cycles.

**Standard**: Store lyophilized standard at -20  $^{\circ}$ C to -80  $^{\circ}$ C for up to 6 months from date of receipt. Aliquot and store the reconstituted standard at -80  $^{\circ}$ C for up to 1 month. Avoid repeated freeze-thaw cycles.

## **GENERAL ELISA PROTOCOL**

#### **Plate Preparation**

1.Dilute the capture antibody to the working concentration in CBS. Immediately coat a 96-well microplate with 100 $\mu$ L per well of the diluted capture antibody. Seal the plate and incubate overnight at 4°C.

2.Aspirate each well and wash with at least 300µl wash buffer, repeating the process two times for a total of three washes. Complete removal of liquid at each step is essential for good performance. After the last wash, remove any remaining wash buffer by inverting the plate and blotting it against clean paper towels. 3.Block plates by adding 300  $\mu$ L of blocking buffer to each well. Incubate at room temperature for a minimum of 1 hour.

4.Repeat the aspiration/wash as in step 2. The plates are now ready for sample addition.

#### Assay Procedure

1.Add 100  $\mu$ L of sample or standards in sample dilution buffer per well. Seal the plate and incubate 2 hours at room temperature.

2.Repeat the aspiration/wash as in step 2 of plate preparation.

3.Add 100 µL of the detection antibody, diluted in antibody dilution buffer, to each well. Seal the plate and incubate 1 hour at room temperature.

4.Repeat the aspiration/wash as in step 2 of plate preparation.

5.Add 200  $\mu$ L of substrate solution to each well. Incubate for 20 minutes at room temperature ( **if substrate solution is not as requested, the incubation time should be optimized** ). Avoid placing the plate in direct light.

6.Add 50  $\mu$ L of stop solution to each well. Gently tap the plate to ensure thorough mixing.

7. Determine the optical density of each well immediately, using a microplate reader set to 450 nm.

## **CALCULATION OF RESULTS**

•Calculate the mean absorbance for each set of duplicate standards, controls and samples. Subtract the mean zero standard absorbance from each.


•Construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph.

•To determine the concentration of the unknowns, find the unknowns' mean absorbance value on the yaxis and draw a horizontal line to the standard curve. At the point of intersection, draw a vertical line to the x-axis and read the concentration. If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

•Alternatively, computer-based curve-fitting statistical software may also be employed to calculate the concentration of the sample.

### **TYPICAL DATA**

This standard curve is only for demonstration purposes. A standard curve should be generated for each assay.



Concentration (pg/ml)	Zero standard subtracted OD
0	0.000
62.5	0.038
125	0.078
250	0.143
500	0.303
1000	0.621
2000	1.274
4000	2.573

#### PERFORMANCE CHARACTERISTIC

#### SENSITIVITY

The minimum detectable dose of Human Contactin 5 / CNTN5 was determined to be approximately 62.5 pg/ml. This is defined as at least three times standard deviations above the mean optical density of 10 replicates of the zero standard.

## TROUBLE SHOOTING

Problems	Possible Sources	Solutions			
	Incorrect or no Detection Antibody was added	Add appropriate Detection Antibody and continue			
No signal	Substrate solution was not added	Add substrate solution and continue			
	Incorrect storage condition	Check if the kit is stored at recommended condition and used before expiration date			
Poor Standard Curve	Standard was incompletely reconstituted or was inappropriately stored	Aliquot reconstituted standard and store at -80 $^\circ\!\!\mathbb{C}$			
	Imprecise / inaccurate pipetting	Check / calibrate pipettes			
	Incubations done at inappropriate temperature, timing or agitation	Follow the general ELISA protocol			
	Background wells were contaminated	Avoid cross contamination by using the sealer appropriately			
	The concentration of antigen in samples was too low	Enriching samples to increase the concentration of antigen			
Poor detection value	Samples were ineffective	Check if the samples are stored at cold environment. Detect samples in timely manner			
High Background	Insufficient washes	Use multichannel pipettes without touching the reagents on the plate			
	insuncient wasnes	Increase cycles of washes and soaking time between washes			
	TMB Substrate Solution was contaminated	TMB Substrate Solution should be clear and colorless prior to addition to wells			
	Materials were contaminated.	Use clean plates, tubes and pipettes tips			
Non-specificity	Samples were contaminated	Avoid cross contamination of samples			
	The concentration of samples was too high	Try higher dilution rate of samples			

ELISA Plate Template	12								
	11								
	10								
	6								
	œ								
	~								
	9								
	Ŋ								
	4								
	m								
	7								
	H								
		۷	8	U	٥	ш	L	U	Ŧ

http://www.sinobiological.com Tel:+86-400-890-9989

Fax :+86-10-58628220

## Human Contactin 5 / CNTN5 ELISA Pair Set Notes