# HIV-1 (group M, CRF07\_BC) GP120 ORF mammalian expression plasmid, N-Flag tag



#### Catalog Number: VG11233-NF

## **General Information**

Gene :	HIV-1 GP120
Official Symbol :	GP120
Synonym :	gp120
Source :	HIV
cDNA Size:	1521bp
RefSeq :	JQ905580.1
Plasmid:	pCMV3-Flag-GP120

### Description

Lot : Please refer to the label on the tube

**Sequence Description :** 

Identical with the Gene Bank Ref. ID sequence.

Restriction site: Kpnl + Xbal (6kb + 1.54kb)

Vector : pCMV3-N-FLAG

Shipping carrier :

Each tube contains approximately 10 µg of lyophilized plasmid.

#### Storage :

The lyophilized plasmid can be stored at ambient temperature for three months.

#### **Quality control :**

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

#### Sequencing primer list :

pCMV3-F:	5' CAGGTGTCCACTCCCAGGTCCAAG 3'	
pcDNA3-R :	5' GGCAACTAGAAGGCACAGTCGAGG 3'	
Or		
Forward T7 :	5' TAATACGACTCACTATAGGG 3'	
ReverseBGH :	5' TAGAAGGCACAGTCGAGG 3'	

pCMV3-F and pcDNA3-R are designed by Sino Biological Inc. Customers can order the primer pair from any oligonucleotide supplier.

#### **Plasmid Resuspension protocol**

1.Centrifuge at  $5,000 \times g$  for 5 min.

2.Carefully open the tube and add 100  $\mu l$  of sterile water to dissolve the DNA.

3.Close the tube and incubate for 10 minutes at room temperature.

4.Briefly vortex the tube and then do a quick spin to concentrate

the liquid at the bottom. Speed is less than  $5000 \times g$ .

5.Store the plasmid at -20 ℃.

#### The plasmid is ready for:

- Restriction enzyme digestion
- PCR amplification
- E. coli transformation
- DNA sequencing

# *E.coli* strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5 $\alpha$  and TOP10F<sup>'</sup>.</sup>

# HIV-1 (group M, CRF07\_BC) GP120 ORF mammalian expression plasmid, N-Flag tag



Catalog Number: VG11233-NF

### **Vector Information**


All of the pCMV vectors are designed for high-level stable and transient expression in mammalian hosts. High-level stable and non-replicative transient expression can be carried out in most mammalian cells. The vectors contain the following elements:

•Human enhanced cytomegalovirus immediate-early (CMV) promoter for high-level expression in a wide range of mammalian cells.

Hygromycin resistance gene for selection of mammalian cell lines.

• A Kozak consensus sequence to enhance mammalian expression.

#### Physical Map of Plasmid :



Vector Name	pCMV3-N-FLAG
Vector Size	6098bp
Vector Type	Mammalian Expression Vector
Expression Method	Constitutive, Stable / Transient
Promoter	CMV
Antibiotic Resistance	Kanamycin
Selection In Mammalian Cells	Hygromycin
Protein Tag	FLAG