

Human Liver cDNA

Catalog #: HD-314

Quantity: 30 Reactions

Storage Conditions: Store at -20°C. It is good for one year of the date of purchase if stored properly.

Applications: The cDNA is primed with oligo dT primer and is ideal for gene expression analysis by PCR, characterization of alternative splicing of mRNA, and Gene cloning and target sequencing.

Quality Control: The PCR-Ready cDNA is functionally tested with the control primers in the recommended PCR reaction.

Description: The PCR-ready first strand cDNA is synthesized from high quality RNA isolated from adult human normal healthy tissues. Total RNA used for cDNA synthesis is isolated by modified guanidine thiocyanate techniques and treated with RNase-free DNase. The total RNA was primed with oligo dT primer and reverse transcribed by a reverse transcriptase enable synthesis of full length cDNA up to 8.9kb. Use μ l cDNA for each PCR reaction.

The amplification conditions used for amplification of beta-actin as a positive control in a volume of 100 μ l:

cDNA	1 μ l
10X polymerase reaction buffer	10 μ l
MgCl ₂ , 25mM (2nM final)	7.8 μ l
Nucleotide Mix, 10mM (0.2mM final)	2.0 μ l
Upstream beta-Actin primer (100 μ M)	1 μ l
Downstream beta-Actin primer (100 μ M)	1 μ l
Taq Polymerase (5 units)	1 μ l
Water	<u>76.2μl</u>
	100 μ l

Program of amplification:

Denaturation 94°C for 2 minutes

25 cycles:

Denaturation 94°C for 1 minute

Annealing 60°C for 1 minute

Extension 72°C for 2 minutes

Final extension 72°C for 5 min

Hold 4°C

Intended use of the products: All products sold by Zyagen are intended for research use only and must not be used for medical diagnostic or drug purposes. Zyagen products may not resell, modified for resale or used to manufacture commercial products without written approval from Zyagen.

Limited product warranty: Zyagen has a total quality assurance policy. If any product does not meet the specification stated in the product data sheet will be replaced at no charge. This warranty limits Zyagen liability to the replacement of the product only. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Zyagen. Zyagen shall not liable for any incidental, indirect, special or consequential damages, even if it is aware of the possibility of such damages.

Zyagen Laboratories, San Diego, California, USA. Telephone: +1 858 546 0720; Fax: +1 858 546 0736; Technical support e-mail address: zinfo@zyagen.com. For additional product and distributor information, see www.zyagen.com